

The skin rejuvenation sensation



Beautiful skin is not a surprise.

It is the visible sign of a good treatment!

Advantages of the treatment with the COLLAGEN Pro Beauty light wave appliance in combination with the Hyaluron Pro Beauty Fluid from Cosmedico:

- Allows skin rejuvenation at home with visible anti-aging effect.
- Improvement of skin elasticity and suppleness.
- Renewal of collagen and elastic fibers.
- Activation of the metabolism of the skin, promotion of youthful skin.
- Reduction of lines (e.g. around eyes and lips, neckline, cheeks).
- Improved well-being, relaxation effect for body and mind.
- No negative side-effects.

► **Switch on and experience the skin rejuvenation sensation first hand!**

PRO BEAUTY
collagen

Stimulation of natural processes

The COLLAGEN Pro Beauty light wave appliance makes the dream of skin rejuvenation come true. Longterm studies have proven that the formation of the body's own collagen is stimulated with red light of a specific wavelength range. The COLLAGEN Pro Beauty tabletop device from Cosmedico emits light in this ideal wavelength range and thus provides exactly the red radiation that stimulates collagen formation and activates skin renewal.

Skin aging factors

Mainly two factors determine the skin aging process: Firstly, the older the skin, the slower the metabolism of the skin. This leads to a decrease in cellular energy in the fibroblast cells that produce collagen and elastin. Secondly, the ability of the body to produce the necessary nutrients for these cells reduces. The result is a slowed production of collagen and elastin, resulting in a skin with less tension and elasticity and thus in the appearance of the first lines.

The gentle waves of red light of the COLLAGEN Pro Beauty table top stimulate the skin to reactivate the formation of skin-firming collagen and elastin. Also, the photo-bio-stimulation in the deeper layers of the skin activates the metabolism in the skin which has a visible smoothing effect on the skin surface.

Before and after

The result of the regular treatment with the COLLAGEN Pro Beauty device is a tense and firm skin with fewer lines and a fresh, youthful appearance.



Regularity pays off

The regular application of red light waves positively rejuvenates the skin complexion. Recommended treatment duration is 12 minutes per application. At the beginning, a bi-weekly treatment is appropriate. Later, it can be extended to up to 3 applications per week. For a visible anti-aging effect, the red light treatment should be pursued for at least 8 - 10 weeks. After a pause of several weeks, the next treatment period can be started. There is no harm at all if you conduct treatments on a regular basis.

The use of Hyaluron Pro Beauty Fluid in combination with the collagen light wave treatment supports the rejuvenation effect significantly and is highly recommended. Apply the fluid prior to each treatment. Close eyes during treatment. Treatment is found to be very relaxing due to the soft warmth of the red light. Application is suitable for all skin types and ages (the skin aging process starts with 25 years), both male and female.

COLLAGEN Pro Beauty light wave appliance



Hyaluron Pro Beauty Fluid (30 ml)



Contact

JW Sales GmbH
Brand Division Cosmedico Medical Systems
Kölner Strasse 8
D-70376 Stuttgart
Tel. +49 (0)711 54 00 4 - 0 | Fax +49 (0)711 54 00 4 - 55
info@cosmedicomed.de | www.cosmedicomed.de